
LiveCodim

Super-resolution for every researcher's desk



From conventional to super-resolution microscopy

LiveCodim is a universal, super-resolution imaging platform, designed to interface with any standard fluorescence microscope.







Telight LiveCodim


LiveCodim is the solution for live imaging with high resolution and low phototoxicity.

LiveCodim is a universal add-on desiged for any fluorescence standard microscope. It offers widefield, confocal, and super-resolution modalities to provide a complete solution, especially suitable for retrofitting.

The technology is based on the patented Conical Diffraction Microscopy - a powerful beam shaper, generating controlled and localized light patterns used to scan biological samples with very low phototoxicity and negligible photobleaching.







Widefield

Confocal

Super-resolution

Cell division

COS7 cell, labeled for actin (purple), mitochondria outer membrane (Green), and tubulin (red). Image shows metaphase using LiveCodim widefield, confocal, or super-resolution (CODIM) modalities.

Key advantages

- High-performance SUPER-RESOLUTION
 MULTICOLOR IMAGING | 4 fluorescence
 channels
- RESOLUTION better or equal 120 nm XY |
 Increased resolution using ligh structuration
 and scanning point approach
- STANDARD SAMPLE PREPARATION |
 Standard fluorescence microscopy sample preparation workflow
- PENETRATION DEPTH | Z-stack ability and "suppression" of signal from other planes

 \bigcirc

SEAMLESS INTEGRATION with all fluorescence microscopes

 \bigcirc

LIVE-CELL IMAGING | Low photo-toxicity and photobleaching enabling in vivo SR study

 \bigcirc

Streamlined and **EASY-TO-USE** application | Autonomous processing algorithm with no need for user-adjustment of many parameters

Specifications Telight Live



Imaging

Modalities

Conical Diffraction SuperResolution microscopy Widefiels epifluorescence Quantitative 5D Live-Imaging

Software features

Automated adaptative SR image processing Fast switching of acquisition mode (Widefield - Confocal - SR CODIM) Guided calibration wizard Intuitive user-friendly interface

Achievable resolution

2D 120nm (Rayleigh) Confocal axial resolution

Field of view

SR CODIM mode: $30 \times 30 \ \mu m$ Confocal mode: $30 \times 30 \ \mu m$ without stitching, $60 \times 60 \ \mu m$ with stitching

Imaging Colors

4 laser line in standard (extra line up to request)

Samples

Standard sample preparation workflow

Hardware

O Dimensions and compatibility

57mm x 47mm x 145mm Inverted Stand from Zeiss, Leica, Nikon and Olympus Universal OME-TIFF File format

Camera

Photometrics Kinetix 95% peak QE 1.2 electrons rms read noise at standard scan

⊘ LED Illumination

CoolLED PE4000 (16 wavelenghts from 365nm to 770 nm)

Laser lines

Violet 405nm Diode Laser (100mW) Blue 488nm Diode Laser (100mW) Green 561 Laser (100mW) Red 640nm Diode Laser (100mW)

Conditions

Objective requirements (for maximum resolution)

60/63x 1.4 (or more) NA (oil immersion) - provided upon request

Right or Left Port C-Mount Compatible Inverted Stand From Zeiss, Leica, Nikon or Olympus Back-Port for LED Illumination Automated XYZ Stage Recommended: Hardware autofocus system (Zeiss DF, Nikon PF, Leica AFC, Olympus IMS) Automted Filter cube Turret

Optical Table