

Luminex[®]

Amnis[®] Imaging Flow Cytometers

Microscopy in Flow




Spanning the Research Disciplines in the Life Sciences.


Microscopy offers detailed cellular images and morphologic information, which are useful scientific tools for the study of cell function. However, the interpretation of microscopy images can be subjective, qualitative, and laborious.

Flow cytometry is excellent for quantitative phenotyping and yields statistically robust results by rapidly interrogating large numbers of cells. However, flow cytometry lacks any ability to image, so sub-cellular localization and functional studies are difficult at best.


By combining the speed, sensitivity, and phenotyping abilities of flow cytometry with the detailed imagery and functional insight of microscopy, the Amnis® ImageStream®X Mk II and Amnis® FlowSight® Instruments overcome the limitations of both techniques and open the door to an extensive range of novel applications.




Immunology




Oncology




Biochemistry




Drug Discovery




Stem Cell Biology




Hematology



Microbiology




Virology



Nanotechnology




Toxicology




Parasitology



Oceanography



Phycology



Small Particle Analysis

Amnis® FlowSight® Imaging Flow Cytometer



- Capable:** Applicable to every research discipline
- Sensitive:** Camera-based detection dramatically increases resolution over traditional flow cytometry
- Affordable:** Smaller footprint with configurations for any lab focus and budget
- Powerful:** Characterizes populations by virtually any visual or fluorescent attribute

Amnis® ImageStream®X Mk II Imaging Flow Cytometer



- High-throughput:** Analyzes thousands of cells per second at up to 60X magnification
- Intuitive:** Simple user interface with real time plotting and gating
- Adaptable:** Can be configured with 1 to 6 lasers
- Boundless:** Variable magnification images small particles and your largest cells

Powerful Flow Cytometry.

The ImageStream[®] Mk II and FlowSight Systems deliver multiple images of every cell in flow, including brightfield, darkfield (SSC), and up to 10 fluorescent markers at high speed. The ImageStream[®] Mk II camera operates with a pixel size of 0.1/0.25/1 μm^2 with 60X/40X/20X magnification, respectively, allowing visualization of fluorescence location from the membrane, cytoplasm, subcellular organelles, and nucleus at high resolution. The FlowSight System operates at 20X magnification with a 1 μm^2 pixel.



The innovative design of Amnis Flow Cytometers increases signal and minimizes noise to provide exceptional photonic sensitivity. Design details like a dedicated side scatter laser, adjustable laser intensities, and brightfield imagery for the direct measurement of cell size allow the systems to resolve cell populations more effectively than far more expensive cytometers. The ease of use, outstanding performance, and imagery of each cell meet the needs of flow cytometry novices and experts alike.

Beyond forward and side scatter

Traditional flow cytometers do an admirable job of using low-resolution scattering characteristics to approximate size and intracellular granularity. Amnis imaging flow cytometers produce familiar ‘size vs. complexity’ scatter plots, but with the power of 20X magnification—or more—can report absolute rather than relative cell size by measuring the actual diameter of objects in brightfield images.

Multichannel immunophenotyping


Immunophenotyping requires multiple fluorescence channels in addition to dual scatter. Below is a 6-color immunophenotype of human PBMC using antibodies against CD3, CD4, CD14, CD16, CD45, and CD123, plus DAPI. The arrangement of detection channels, available laser options, and automated compensation wizard allows for the straightforward separation of complex cell populations.



Sensitive and Flexible for Any Research Need.

Exceptional fluorescence sensitivity


The patented architecture of Amnis imaging flow cytometers provides extraordinary fluorescence sensitivity across the visible spectrum, outperforming other imaging devices. The four plots below demonstrate the ability of the FlowSight Instrument to discriminate all intensities in the Spherotech 8-peak calibration bead set, across the spectrum from FITC to PE-Cy7. Note the distinct peak separation, low coefficients of variance (CVs), and high sensitivity from the FITC to the PE-Cy7 channels.



FlowSight[®] 12 channel imagery of 3-micron diameter Spherotech 8-peak Rainbow beads.

5-Part white blood cell differential

Because of its exceptional sensitivity, the FlowSight System excels at the resolution of mixed sub-populations in heterogeneous samples. Human peripheral blood mononuclear cells (PBMC) are partitioned into 5 distinct populations using CD45 expression and side scatter intensity. High fluorescence sensitivity and tight coefficients of variance (CVs) resolve monocytes (green) from lymphocytes (blue) and facilitate the detection of rare basophils (white). The dedicated side scatter laser clearly resolves eosinophils (yellow) from neutrophils (orange).




Sensitive and Flexible for Any Research Need.

Images of every cell


The FlowSight and ImageStream[®] Mk II Instruments operate like conventional flow cytometers, but also provide imagery of every cell. Powerful and intuitive analysis software seamlessly links quantitative data to images:

- Click on a dot in any plot to see its corresponding image
- Click on a bin in any histogram to view every cell in that bin
- Draw gates on dot plots and view the resulting populations to validate results




With imaging capabilities, you'll never wonder about outliers or whether your gates are in the right place, as shown in the example above. Once you've drawn a gate on a plot, you can click inside and out to determine if it's in the right place, as shown in the example to the right. With visual feedback, you can optimize gate size, shape, and position for better data quality.

Granulocytes




Monocytes



Double Positives



Doublet Artifacts



Data Acquisition Software.

Fast and easy.

INSPIRE[™] Software offers powerful, image-based gating and real time fluorescence compensation.

- 1. Instant Population Viewer:** Every population is added to a pull-down list as soon as you draw a gate. Simply select a population of interest from the list to view the corresponding cells during data acquisition.
- 2. Image Gallery:** Imagery of cells of interest appear in the gallery as they are acquired, allowing you to inspect morphology, assess staining patterns, and optimize laser power settings.
- 3. Instrument Status at a Glance:** Convenient gauges, indicators, and text alerts provide continuously updated instrument operational status.
- 4. Real Time Intensity Compensation:** An easy to use compensation wizard quickly guides you through the setup of multi-color compensation matrices.
- 5. Gating Without Guesswork:** Gates are easily drawn using graphical tools, and verified for accuracy by visual inspection of gated cells.
- 6. Efficient Sample Handling:** Up to 95% of the sample volume is utilized, facilitating the analysis of rare cells. Unused sample can be recovered for further analysis.
- 7. Intuitive Acquisition:** A simple and intuitive user interface provides complete control of sample acquisition settings and data storage criteria.
- 8. Familiar Dot Plots and Histograms:** Data plots are updated in real time, just as with conventional flow cytometers. Unlike conventional cytometers, you can also plot morphologic parameters such as Area, Cell Width, Cell Height, Aspect Ratio, and others.

INSPIRE[™] Software



Software That Turns Data Into Understanding.

IDEAS® Software combines image analysis, statistical rigor, and visual confirmation in an easy to use package

- 1. Inspect Your Populations:** The Image Gallery allows you to see every image of every cell or perform a “virtual cell sort” to inspect and validate the cells within a specific population.
- 2. Images for Every Dot:** Every dot in every scatter plot is linked to the corresponding cell imagery. Simply click on a dot to see the associated cell images or vice-versa.
- 3. Graphical Population Definitions:** Define populations using familiar graphical tools and combine them with logical functions.
- 4. Comprehensive Population Statistics:** Characterize your cell populations with a wide range of statistical metrics to reveal differences in cell morphology, phenotype, and function.
- 5. Flexible Image Display Tools:** Create composite images, pseudo-color representations, and a host of other image transformations for reporting and publication.
- 6. Graph What You See:** Virtually anything you see in the imagery can be plotted as a histogram or dot plot. Hundreds of parameters are calculated for every cell, including fluorescence intensity, fluorescence location, cell shape, cell texture, and numerous other morphologic and photometric features.

IDEAS® Software



A Wealth of Applications.


Any application you can imagine.

Featured applications


The applications detailed on the following pages demonstrate the types of studies that can be performed using the ImageStream[®] Mk II and FlowSight Instruments with their powerful companion IDEAS image analysis software.

Any application you can imagine


The ImageStream[®] Mk II and FlowSight Systems are designed to be general-purpose platforms for cellular studies and are not limited to the applications illustrated in this brochure.




Cell Signaling




Cell Cycle and Mitosis




Internalization and Co-localization




Surface and Intracellular Co-localization




DNA Damage and Repair



Cell-cell Interaction




Shape Change and Chemotaxis




Immunological Synapse



Micronucleus Counting




Oceanography



Stem Cell Biology



Parasitology




Microbiology


Quantifying Nuclear Translocation...

20X resolution tells the story


Translocation of NFkB from the cytoplasm to the nucleus of the cell is a key event in the response to the presence of cell stressors. Only imaging flow cytometers can analyze translocation quantitatively, in thousands of cells. For this data, the 20X objective of the FlowSight System is used to locate NFkB in relation to 7-AAD fluorescence from the nucleus in untreated THP-1 cells and cells stimulated with lipopolysaccharide (LPS). The similarity feature of the IDEAS Software produces a score for every cell, quantifying the co-localization of NFkB and 7-AAD.



THP-1 Control (no LPS)
Mean similarity score = 0.4




THP-1 + 1 μg/mL LPS
Mean similarity score = 3.2




...with Flourescent Image Similarity.

A closer look at NFkB signaling with 60X magnification


Here, THP-1 cells stimulated or not with LPS and stained with anti-NFkB and 7-AAD to counterstain the nucleus were collected on the ImageStream[×] Mk II System using the 60X objective. The IDEAS Software similarity feature demonstrates binning of samples consistent with the FlowSight histogram, and establishes the quality of visual detail that the ImageStream[×] Mk II System can provide for studies when greater details are needed.



THP-1 Control (no LPS)
Mean similarity score = 0.2



THP-1 Control (no LPS)
Mean similarity score = 3.8



Amnis® Spectral Imaging Channels And Corresponding Fluorophores

Laser	Fluorophore	Ex	Em	☀	Fluorophore	Ex	Em	☀	Fluorophore	Ex	Em	☀	Fluorophore	Ex	Em	☀	Fluorophore	Ex	Em	☀	Fluorophore	Ex	Em	☀	FlowSight [®]	ImageStream [®] Camera 1	ImageStream [®] Camera 2	
375 (with installed 405)	CH 1				CH 2				CH 3				CH 4				CH 5				CH 6				●	●	●	
	Ch1/Ch9 BF *or*				QD525				eFluor565 NC				eFluor625 NC				eFluor700 NC				QD800							
	Alexa Fluor® 350				350-450				UV - 405				UV - 405				UV - 405				350-450							
	BV421™				525				565				625				705				800							
	Cascade Blue				5				585				625				705				5							
	DAPI																											
488	Pacific Blue																								●	●	●	
					Alexa Fluor® 488				Cy3				AldeRed				7-AAD				PE-Alexa Fluor® 750							
					BODIPY FI				514				488				546				496,565							
					DiO				557				PE-Alexa Fluor® 610				DRAQ5				496,565							
					DyLight™ 488				592				PE-Texas Red®(ECD)				FuraRed-lo				775							
					FITC				578				RFP				LDS751				774							
					494				584				555				543				712							
					520								PE-Alexa Fluor® 647				496,565											
					509								PE-Cy5				496,565											
					LysoTracker Green								PE-Cy5.5				496,565											
					MitoTracker Green								PerCP				482											
					PKH2 & PKH67								PerCP-Cy5.5				482											
					504								PI				535											
					Rhodamine 110																							
					SYBR® Green																							
561					Syto13 (DNA/RNA)																				●	●	●	
					YFP																							
									Alexa Fluor® 546				Alexa Fluor® 568				7-AAD				PE-Alexa Fluor®750							
									CellMask/Tracker				DyLight™594				DRAQ5				496,565							
									522				593				646				496,565							
									549				618				712				774							
									565				613				712											
									DSRed				PE-Alexa Fluor®610				PE-Alexa Fluor® 647				669							
785																									●	●	●	
Ch width	435-480				480-560				560-595				595-642				642-745				745-780				Ch width			
Bandpass*	(457/45)				(528/65)				(577/35)				(610/30)				(702/85)				(762/35)				Bandpass*			
375 (with 405 not installed)	CH 7				CH 8				CH 9				CH 10				CH 11				CH 12				●	●	●	
	Alexa Fluor® 350				eFluor525 NC																							
	BV421™				UV - 405																							
	Cascade Blue				350-450																							
	DAPI				525																							
	Hoechst				5																							
405	Pacific Blue																								●	●	●	
					Alexa Fluor® 430								eFluor625 NC				eFluor650 NC				QD800							
					434								UV- 405				UV- 405				350-450							
					BV510™								625				650				800							
					405								625				705				5							
					Cascade Yellow																							
					Pacific Orange																							
					410																							
592					Pacific Orange																				●	●	●	
					410																							
					551																							
					DAPI																							
					410																							
					DyLight™405																							
					420																							
					Hoescht																							
642					LIVE/DEAD Violet																				●	●	●	
					416																							
					451																							
					Pac Blue																							
					455																							
785																									●	●	●	
Ch width	435-505				505-570				570-595				595-642				642-745				745-780							Ch width
Bandpass*	(457/45)				(537/65)				(582/25)				(610/30)				(702/85)				(762/35)				Bandpass*			

Quantitative Imaging and Robust Population Statistics.

Quantitative imaging means Luminex's imaging flow cytometers have a powerful and intuitive image processing package with thousands of analysis parameters and optimized analysis wizards for many common image-based applications, including nuclear translocation, shape change, internalization, and apoptosis.

Objective, quantitative image analysis on large numbers of cells is backed by a large set of statistical parameters for data reporting.

Select the wizard to use for analysis:

Open File

Display Properties

Begin Analysis

Feature Finder

Apoptosis

Cell Cycle - Mitosis


Co-localization

Internalization

Nuclear Localization

Shape Change

Spot




File	Count All	Count Focus	Count Singles	Count Positive	Mean Similarity	Std Dev Similarity
TNFa_0ng_2_2016.daf	10000	4903	4265	3740	0.34	0.71
TNFa_0-05ng_3_2016.daf	10000	4621	4060	3635	1.28	0.81
TNFa_0-1ng_4_2016.daf	10000	4280	3739	3365	1.90	0.82
TNFa_0-5ng_5_2016.daf	10000	4861	4167	3516	2.68	0.66
TNFa_1ng_6_2016.daf	10000	3811	3311	2910	2.72	0.63
TNFa_2-5ng_7_2016.daf	10000	3893	3425	3070	2.80	0.58
TNFa_5ng_8_2016.daf	10000	4162	3685	3180	2.80	0.52
TNFa_7-5ng_9_2016.daf	10000	4361	3782	3387	2.82	0.58
TNFa_10ng_10_2016.daf	10000	4005	3456	2988	2.90	0.55



Internalization Identifies Trogocytosis.

20X objective for a wider field of view

The FlowSight Instrument is optimized for imaging large objects such as epithelial cells, macrophages, neutrophils, fibroblasts, and even large eukaryotic parasites. Here, *Entamoeba histolytica* demonstrates amoebic trogocytosis of immune cells. Following attachment to Jurkat cells, the FlowSight Instrument measures every *E. histolytica* expressing Jurkat markers internalized or on their surface.




Data courtesy of Dr. Katherine Ralston, UC Davis.



Co-localization and Trafficking.

The ImageStream[®] Mk II Instrument greatly improves co-localization studies by combining the rapid collection of large numbers of cell images with objective measurement of the similarity of bright image details.





Lysosomal trafficking of CpGB within pDC is quantified using the Internalization (Y-axis) and the Bright Detail Similarity (X-axis) scores. Representative merged images of pDC (orange), CpGB (red), and lysosomes (green) are shown at 40X magnification cells within the lower left region of the plot and have surface-bound CpGB. As CpGB molecules enter the pDC, the Internalization score increases (upper left region). Once the CpGB traffics to the lysosomes, the similarity between the CpGB and lysosome image pair increases (upper right region).

Data courtesy of Dr. Patricia Fitzgerald-Bocarsly, University of Medicine and Dentistry, New Jersey.

Apoptosis and Necrosis.


Apoptosis and necrosis detection by image analysis

The apoptosis wizard analyzes the nuclear morphology and brightfield image contrast of each cell to detect apoptosis in any sample containing a nuclear stain.




Necrosis versus apoptosis

Conventional flow cytometers can use membrane-impermeant dyes to identify dead or dying cells that have lost membrane integrity. However, it can be difficult to determine if cell death is via apoptosis or necrosis. The FlowSight System simplifies this determination by revealing the nuclear morphology of every cell. As shown in this sample of THP-1 cells labeled with propidium iodide, the nuclei of necrotic cells have normal nuclear morphology, while the nuclei of apoptotic cells are shrunken and fragmented.



Autophagy.


During autophagy, cytoplasmic LC3 is processed and recruited to the outer membrane of autophagosomes. Cells undergoing autophagy can be identified by visualizing LC3 puncta and enumerating the spots within each cell using the Spot Count feature of the IDEAS Software package.



The IDEAS image processing software included with the ImageStream[®] Mk II System determines the Spot Count of every cell. In this example, cells with varying number of LC3-RFP (red) spots are shown with their corresponding Spot Count.

Example: Autophagy in the human osteosarcoma cell line U2OS

A serum starvation step can induce autophagy in U2OS cells. This data demonstrates how the IDEAS feature Spot Count can be used to quantify the textural variation between healthy cells and those undergoing autophagy.




U2OS RFP-LC3 human osteosarcoma reporter cell line was starved for 4 hours at 37°C. Both the control and starved samples were supplemented with a degradation inhibitor. FlowCollect RFP-LC3 Reporter Autophagy Kit (Catalog No. FCCH100183).

Morphology.


Change in cell shape is correlated with change in function, particularly in the case of macrophage activation, stem cell differentiation, and cellular response to drugs. The ImageStream[®] Mk II Instrument measures cell shape using powerful, pre-defined features in the IDEAS image analysis software. One such feature is the Circularity score.

The Circularity score is a measure of how much the cell radius varies. Round cells (left) have high Circularity scores while irregularly shaped cells (right) have low Circularity scores.



Example: Shape change in primary monocytes

Chemoattractant MCP-1 induces monocyte shape change and migration to sites of inflammation, as evidenced by the significant decrease in the Circularity score of the MCP-1 treated sample relative to the untreated control. In contrast, treatments that reduce inflammatory response – such as drugs for autoimmune disorders – result in an increase in Circularity scores.




Microalgae.

Quintessential Cell Interactions at the Immunological Synapse.


Mixed cultures of microalgae

The images below demonstrate microalgae identification in mixed cultures using morphological parameters and the ImageStream[®] Mk II Instrument at 40X magnification.




Microalgae quality control


The images below demonstrate detection of bacterial contamination, cellular debris, and clusters in mixed culture of microalgae. A mixed culture of *T. pseudonana* and *C. cryptica* contaminated with bacteria was analyzed on the ImageStream[®] Mk II System at 60X magnification.



FlowSight 20X images




SEB Dose Response Curve




Raji B cells were exposed to SEB (0-20 μg/mL) and incubated with human primary T cells.

Take the Analysis Even Further with Higher Resolution.


- T:APC conjugates are easily identified using morphological features
- The point of cell-cell contact is identified using a mask (cyan overlay)
- Actin accumulation within the mask confirms formation of an immunological synapse
- All T cells are then identified either in conjugates or not
- NFkB translocation is measured in the T cells specifically



ImageStream[×] Mk II 60X images




Modular Options for the **FlowSight** and **ImageStream[×] Mk II** Instruments.




Additional excitation lasers

The 488 nm blue laser comes standard with the FlowSight and ImageStream[×] Mk II Instruments. Adding excitation lasers increase experimental flexibility by permitting a broader palette of fluorescent markers. All lasers are intensity adjustable to ease protocol development.




12 channels of detection

Up to 12 high resolution image channels are available with the addition of an optional second camera and associated optics for the ImageStream[×] Mk II System. Twelve channels are standard on the FlowSight Instrument.




Multi-well plate AutoSampler

The AutoSampler option enhances productivity with unattended sample loading from 96-well plates. The fully integrated AutoSampler option greatly facilitates dose response and time course studies.



MultiMag

The MultiMag option for the ImageStream[×] Mk II System provides 60X and 20X objectives on a motorized stage, in addition to the standard 40X objective. The 60X objective offers greater resolution for the morphologic analysis of cells as small as yeast and bacteria, while the 20X objective offers a 120 micron wide field of view for very large cells.



EDF: Extended depth of field

The EDF option incorporates Wavefront Coding technology from OmniVision CDM Optics, which is a combination of specialized optics and unique image processing algorithms, to project all structures within the cell into one crisp plane of focus. Ideal for automated FISH spot counting.

Option	FlowSight [®]	ImageStream [×] Mk II
Additional Excitation Lasers	Standard 488 Option 405, 561, 642	Standard 488 Option High Power 488, 375, 405, 561, 592, and 642
12 Channels of Detection	Standard	6 Standard High Resolution 12 Channel Option
Multi-well Plate Autosampler	96-Well Plate	96-Well Plate
MultiMag	Not Available	40X Standard; 20X and 60X Option
EDF: Extended Depth of Field	Not Available	Available

Progressive Engineering...

FlowSight Instrument specifications

Performance Characteristics	Magnification 20X
Numeric aperture	0.6
Pixel size	1.0 x 1.0 µm
Field of view	60 x 256 µm
Imaging rate	4,000 cells/sec

Sample characteristics

- **Volume** - 20-200 µL
- **Utilization Efficiency** - Up to 95% of sample

Automated instrument operations

- Start up and shut down
- Sample load and acquisition
- Laser alignment, focus adjustment, calibration, and self-test

Operational requirements

- 400W, 100-240 VAC, 50/60 Hz
- No external air or water necessary

Physical characteristics

- 18" W x 18" H x 25" D in (457 mm x 465 mm x 635 mm)
- 135 lbs. (61 kg)

Illumination

- **Excitation** - Standard: 488 nm; Optional: 405 nm, 561 nm, and 642 nm
- **Side scatter** - 785 nm standard
- **Brightfield** - Multi-channel



...Advances Performance.

ImageStream[×] Mk II Instrument specifications

Performance Characteristics	Magnification		
	40X	60X	20X
Numeric aperture	0.75	0.9	0.5
Pixel size	0.5 x 0.5 µm	0.3 x 0.3 µm	1.0 x 1.0 µm
Field of view	60 x 128 µm	40 x 170 µm	120 x 256 µm
Imaging rate	2,000 cells/sec	1,200 cells/sec	5,000 cells/sec

Sample characteristics

- **Volume** - 20-200 µL
- **Utilization Efficiency** - Up to 95% of sample

Automated instrument operations

- Start up and shut down
- Sample load and acquisition
- Laser alignment, focus adjustment, calibration, and self-test

Operational requirements

- 450W, 100-240 VAC, 50/60 Hz
- No external air or water necessary

Physical characteristics

- 35" W x 26" H x 25" D in (889 mm x 660 mm x 635 mm)
- 400 lbs. (182 kg)


Illumination

- **Excitation** - Standard: 488 nm; Optional: High Power 488, 375 nm, 405 nm, 561 nm, 592 nm, and 642 nm
- **Side scatter** - 785 nm standard
- **Brightfield** - Multi-channel



The Path to Scientific Enlightenment...

...passes through the Amnis® multi-spectral decomposition element, which enables simultaneous collection of brightfield, laser scatter, and multiple fluorescent images per cell.



Ordering Information

Product Name	Part Number
Instruments	
Amnis® FlowSight® Flow Cytometer	100370
Amnis® ImageStream® Mk II Flow Cytometer	100220
Reagents	
Amnis® SpeedBead® Kit	400041
FlowSight® Calibration Beads	400300
Kits	
Amnis® NFkB Translocation Kit	ACS10000
Amnis® Protein Aggregate and Silicone Oil Detection Kit	APH10001
Amnis® Intracellular Staining Kit	ACS10002
Service Plans	
ImageStream® Mk II Intermediate Maintenance and Service Agreement	SLA-ISXMKII-Intermediate
ImageStream® Mk II Complete Maintenance and Service Agreement	SLA-ISXMKII-COMPLETE
ImageStream® Mk II Basic Maintenance and Service Agreement	SLA-ISXMKII-BASIC
FlowSight® Intermediate Maintenance and Service Agreement	SLA-FS-Intermediate
FlowSight® Complete Maintenance and Service Agreement	SLA-FS-COMPLETE
FlowSight® Basic Maintenance and Service Agreement	SLA-FS-BASIC
Amnis® ImageStream® Mk II System IQ/OQ	603222
Amnis® FlowSight® IQ/OQ	603224
Training Options	
ImageStream® Mk II training at Luminex, 3 days - per person	500200
Onsite ImageStream® Mk II or FlowSight® training - FAS 1 day; Up to 5 people	500200-1
Onsite ImageStream® Mk II or FlowSight® training - FAS 2 consecutive days; Up to 5 people	500200-2
Onsite ImageStream® Mk II or FlowSight® training - FAS 3 consecutive days; Up to 5 people	500200-3
Onsite ImageStream® Mk II or FlowSight® training - FAS 4 consecutive days; Up to 5 people	500200-4
Onsite ImageStream® Mk II or FlowSight® training - FAS 5 consecutive days; Up to 5 people	500200-5
Software	
Amnis® IDEAS® Image Analysis Software – 21CFR- enabled	IFC300202
Amnis® INSPIRE™ Software & Amnis® IDEAS® Software – 21CFR-enabled, Amnis® FlowSight®	IFC300203
Amnis® INSPIRE™ Software & Amnis® IDEAS® Software – 21CFR-enabled, Amnis® ImageStream® Mk II	IFC300204
Amnis® IDEAS® 6.3 Image Analysis Software – Single seat license	CN-SW69-01
Amnis® IDEAS® 6.3 Image Analysis Software – 12 seat group license	CN-SW69-12
Amnis® IDEAS® 6.3 Image Analysis Software – Institutional license	CN-SW69-20
Machine Learning Module	CN-SW45-01
Amnis® AI: Computer-aided Image Analysis Software	CN-SW70-01



For more information, please visit luminexcorp.com/flowsight-imaging and luminexcorp.com/imagestreamx-mark-ii.

For Research Use Only. Not for use in diagnostic procedures.

©2020 Luminex Corporation. All rights reserved. Amnis, FlowSight, ImageStream, IDEAS, and SpeedBead are trademarks of Luminex Corporation registered in the U.S. and other countries. INSPIRE is a trademark of Luminex Corporation. All other trademarks are the property of their respective owners.

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europe@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.8036.9888

info@luminexcorp.com

JAPAN

+81.3.5545.7440

info@luminexcorp.com

BR243222